

Description

TAM III micro-calorimeter is an ultra-sensitive heat flow measurement device for real-time monitoring and dynamic analysis of chemical/physical processes, over a period of hours or days. It can be used to determine the onset, rate, extent and energetics of various processes for specimens in small ampoules (3–20 ml) at a constant set temperature (ca. 15°C–150°C) or temperature ramp (0-2K/hr).

Before starting:

- 1) Ensure that the work area is clear of dust, debris, material samples, and any lab components not directly used for the experiment.
- 2) Put on appropriate personal protective gear (safety glasses, lab coat and gloves)

Loading TAM III ampoules using glovebox in Lab 221 for sensitive materials:

- 1) Place approximately 30-60 mg of wet sample (under hexane) in 50 ml glass vial.
- 2) Place glass vial, 20ml disposable crimp seal TAM ampoules, and ampoule cover inside glove box under argon.
- 3) Using balance located inside glove box, accurately record stabilized weight of sample as hexane evaporates.
- 4) Add heavy paraffin oil, CAS: 8042-47-5 to completely immerse the sample.
- 5) Seal TAM ampoules using the crimping tools and attach hanger at the top of ampoule.
- 6) Remove the ampoule from the glove box and carry it to Lab 107 where TAM III is located.

Loading TAM III:

- 1) Select appropriate experimental mode on TAM III interface: isothermal/scanning
Note: Only two samples can be loaded at a given time.
- 2) Input sample weight and other sample details into device software.
- 3) For isothermal mode: select experimental furnace temperature: 30-140°C.
 - a. Allow device to acquire initial baseline heat flow signals for at least 30 min.
 - b. Remove ampoule lifter from device and connect with hanger attached to the sealed ampoule.
 - c. Lower sample into equilibration position using the indicated mark on the ampoule lifter and held it there for 15 min.
 - d. Lower sample into the measuring position and start recording data.
 - e. Stop after desired experimental time (hours/days).
 - f. Remove samples and allow device to acquire final baseline.
- 4) For scanning mode: select desired heating rate: 0-2 K/hr
 - a. Input desired temperature program.
 - b. Remove ampoule lifter from device and connect with hanger attached to sealed ampoule.
 - c. Lower sample directly into measuring position without automatic baseline.
 - d. Allow heat flow signals to stabilize for 3-6 hours.
 - e. Start temperature program.
 - f. Remove samples upon completion of temperature program.

Unloading TAM III ampoules:

- 1) Safely dispose remaining material from ampoule into waste container for reactive powder ensuring appropriate (1:3) dilution with sand.