Wire Ignition – Standard Operating Procedure

Modified 08/14/2013

- a. Turn on the IR pyrometer monitor; note the stabilization time required: 15 min
- 2. Open the software Redlake MiDAS on the computer attached to the High Speed Camera
 - a. Select "Open Camera" in the software
- 3. Straighten wire and cut wire to appropriate length, about 6~7cm.
- 4. Open the wire ignition chamber, load the wire and coat the powder.
 - a. Loosen the holders (both end),
 - b. Place the wire and secure them in the groves etched in the holders
 - c. Tighten up the holder on the right hand side first
 - d. Make sure the wire does not bend.
 - e. Tighten the right holder
 - f. Put a spacer-template between the two holders, and push the left holder against the spacer, then tighten the left holder
 - g. Mix the powder with hexane, put the spacer between the holders again, and use a brush to make a thin coating of powder in the space provided in the spacer-template.
 - h. Lock the chamber
- 5. Select Record on the MiDAS software
 - a. Set the Frame Rate to 25 fps
- 6. Adjusting camera:
 - a. Set aperture initially at 2.5,
 - b. Move the red light into the frame of the camera,
 - c. Choose zoom ring to zoom in/out to focus on the wire and powder,
 - d. Choose focus ring to focus on the powder. Then change the aperture to 8-11
- 7. Open the file Wire Ignition Evin 1 on computer connected to the pyro-meter.
 - a. Choose signal 0 and channel 1, then click on $\sqrt{}$ button to change signal
 - 1 and channel to 3, and then click $\sqrt{\text{again to save}}$,
 - b. Change the number of scans of acquire to 3000
- 8. Turn in the resistance box according to what current we need.
- 9. After turning lights off, change the Frame Rate to 500 fps.
- 10. Once the red light of the thermal monitor is on,
 - a. Click the \rightarrow button on pyro-meter computer,
 - b. Press Record on the main computer
- 11. Once the green light is on, close the switch. When wire is burned, open the circuit / switch.
- 12. Stop video recording.
- 13. Record the frame rate and temperature when the small light is on. Select a time window in the recorded video that captured the ignition event.
- 14. Save the recorded data.
- 15. Save the data from the pyrometer.
- 16. Once when all the runs are completed, turn off the IR pyrometer monitor.