


c. Press button **LIFT ON/OFF** (1b) to insert your sample into the magnet

2) Create a new data set by typing **new** <<enter> at the command line or click on at the toolbar (or type **edc** <enter> at the command line to start from an existent template)



NAME -- (filename for your compound)
 EXPNO -- suggestion: Set 1 for 1H, 2 for a 13C
 PROCNO -- suggestion: set to 1
 DIR -- set to /BRUKER/TOPSPIN/
 USER -- should be <Your_Professor name>
 SOLVENT -- your deuterated solvent from the pull-down menu
 EXPERIMENT -- "Use current params"
 TITLE -- the title of your sample (can be as descriptive as needed)

3) Spin the sample at 25 Hz: press the button shown in arrow 1c in Figure 1.

4) Read a good shimfile with **rsh** <enter>: Choose the **stdshims_CP-300** file from the list, which is the standard shimfile name for the NMR spectrometers.

5) Type **lock** <enter> at TopSpin command line and choose the appropriate solvent from the lock table.

6) After the lock process is finished, shim **z1**, **z2**, **z3**, and **z4** pressing the buttons shown in Figure 1, arrow 1d.

7) Type **rga** <enter> at the TopSpin command line.

a. Adjust the number of scans - NS - as appropriate if necessary; default for 1H is 16, for 13C is 1024)

- 8) Type **zg** <enter> at the command line to start the experiment.
- 9) Your NMR data is automatically saved **at the end** of the experiment
- 10) Finishing the experiment, press lift button -1a- to eject your sample and insert the “dummy” sample, pressing the lift button again, 1b. Cover the probe with the lid and lock on the “dummy” sample (D₂O).

DATA PROCESSING

It is highly recommended to carry out NMR data processing *off-line*


****Bruker is now offering Topspin software free for Academia: the software can be downloaded at:**

<https://www.bruker.com/service/support-upgrades/software-downloads/nmr/free-topspin-processing/free-topspin-download.html>


Optionally, you may use the MNova software (NJIT's NMR laboratory will have 5 licenses) -email me to obtain your license and instructions for downloading and activation


Processing using TopSpin


- Type **ef; apk** <enter> at the command line to perform Apodization (e), Fourier (f) Transformation, and automatic phase adjustment (apk).
- If needed, further optimize the phase of the spectrum manually, **using the phase correction icon.**


- Click on the phase icon: 


Left-click **0**, hold the mouse button to correct phase **at** the vertical cursor. Left-click **1**, hold the mouse button to correct phase **away** from the cursor. Click on the SAVE icon to exit the phasing window.




- Calibrate your spectrum using a known chemical shift (residual solvent peak or TMS). Click on the icon  at the toolbar. Set the reference and return.

- Type **abs n** to flatten the baseline. There is a manual baseline flattening: 

- Integrate the spectrum manually with  -- This takes you into a subroutine: Integrate the spectrum, save the integrals before exit the window.



- Peak picking can be done with the icon  and drawing a green box around the peaks. Save the peaks before exiting the window (using the same icon as the phasing and integration).



- Save your data on the flash drive (alternatively, zip the NMR data directory and email to you as attachment).