



# New Agilent 7900 ICP-MS Product Information

Redefining ICP-MS.  
*Again!*

# 26 Years of Agilent ICP-MS Innovation

*A New Dimension in ICP-MS Capabilities*

**7900 ICP-MS**



*A new era in ICP-MS analysis*

**8800 ICP-QQQ**



*Bringing ease of use and reliability in ICP-MS*

**7700 Series**



*Setting new levels to ICP-MS Performance*

**7500 Series**



*Enabling routine robust ICP-MS analysis*

**4500 Series**



*Pioneering ICP-MS technology*

**PMS 100**



World's first benchtop system.  
Hyperbolic profile quad, motorized torch XYZ, cool plasma

**#1 selling ICP-MS**

**#1 selling ICP-MS**

HMI, ISIS-DS, 3<sup>rd</sup> Generation ORS  
MassHunter SW

**#1 selling ICP-MS**

World's first ICP-QQQ for unparalleled interference removal and sensitivity

With 10x higher matrix tolerance, 10x wider dynamic range and 10x better single to noise the new Agilent 7900 rewrites the rules on ICP-MS

1987

1994

2000

2009

2012

2014



Agilent Technologies

# Unique Performance of the 7700

## Better matrix tolerance than any other ICP-MS

- Higher plasma temperature (lower CeO/Ce ratio) under standard conditions than any other system

## Best performance with Helium cell gas – eliminates need for reaction gases in all common applications

- ◆ 7700 ORS<sup>3</sup> improvements - removes all polyatomics in He mode, giving accurate results in complex or variable sample types – impossible on ICP-MS systems that use reactive cell gases or mixtures




## Wider dynamic range than any other quadrupole ICP-MS

- ◆ Full **9 orders** dynamic range at the detector – linear to 500ppm without changing conditions or hardware




# How to Improve on the Best ICP-MS?

How about a NEW Interface, Ion Lens & Vacuum system for up to 15x **Signal to Noise** improvement?



How about NEW **detector** technology for unprecedented 11 orders measurement range and Fast TRA Acquisition at 0.1msec dwell time.?



How about a NEW **Ultra HMI** for analysis of up to 25% NaCl brine solutions without dilutions?

How about a NEW **ORS<sup>4</sup>** for ultra fast gas switching in less than 3 seconds?

# Presenting Agilent's Game-Changing 7900 ICP-MS



**We took the world's best-selling, highest performing quadrupole ICP-MS, and made it 10x better!**

# Agilent 7900 ICP-MS Key Performance Highlights

## Rewriting the rules on ICP-MS

**10x better signal to noise (S/N)** – even lower detection limits

- 2x higher sensitivity, 3-5x lower background

**10x higher matrix tolerance** – handles even tougher samples

- HMI is still unique to Agilent. On the new ICP-MS, the optional **Ultra HMI (UHMI)** extends capability to matrix levels of **up to 25%**

**10x wider dynamic range** – increases upper measurement limit

- 7700's 9 orders was best in class and new ICP-MS extends this by at least a further order of magnitude (**up to 11 orders detector linear dynamic range**), allowing % levels to be quantified - a first for ICP-MS

**30x faster detector** –faster transient signal measurement (TRA)

- **0.1ms** integration time means improved single nanoparticle analysis

**Improved productivity**

- New ultra fast **ORS<sup>4</sup>** with less than 3 seconds switching time between modes
- New **ISIS 3** for fast unattended start-up, autotune and sample delivery

# Higher sensitivity & lower background than 7700x

## Comparison to 7700x actual performance

~2x higher sensitivity at mid- and high-mass (~same at low-mass)

	Plasma mode	Sensitivity (7900/7700x)		
		7Li	89Y	205TI
No Gas	Normal (spec. tune)	x1.2	x1.9	x1.9
	Low Matrix	x1.0	x1.6	x1.6
	General Purpose	x1.0	x1.5	x1.8
	HMI-high	x1.3	x1.6	x1.8

3x to 5x lower background

Model	Sensitivity [kcps/ppb]			Off Mass Bkgnd [cps]			S/N [x1E3]		
	7Li	89Y	205TI	7amu	89amu	205amu	7Li	89Y	205TI
7700x	100	320	320	0.7	1.5	5.4	140	210	60
<b>7900</b>	<b>120</b>	<b>620</b>	<b>600</b>	<b>0.2</b>	<b>0.2</b>	<b>1.2</b>	<b>600</b>	<b>3100</b>	<b>500</b>
7900/7700x	1.2	1.9	1.9	0.28	0.13	0.22	4	15	8

Average ~10x better signal to noise (S/N)

# Agilent 7900 Performance Specifications


Numbers in red are better spec. than 7700.

		7700	7900
		Spec	Spec
Sensitivity (Mcps/ppm)	7Li	50	55
	89Y	160	320
	205Tl	80	250
Background (cps)	No gas (9amu)	<2	<1
CeO/Ce		<1.5%	<1.5%
CeO/Ce (HMI-25)		-	<0.5%
Ce <sup>++</sup> /Ce		<3%	<3%
DL (ppt) (No gas)	Be	0.5	0.2
	In	0.1	0.05
	Bi	0.1	0.08
DL (ppt) He mode	As	20	20
	Se	40	40
Short-term stability (20 min)		<3%	<2%
Long-term stability (2 hrs)		<4%	<3%



# Agilent 7900 Capable of Very High Sensitivity

## Tuned like competitor's ICP-MS – CeO/Ce <2.5%



**Uranium** calibration in No Gas mode:

Ultra-high sensitivity  
**1.38 GHz/ppm**



Ultra-low background:  
**MDL: 1.3ppq; BEC: 0.48ppq;**

Calibration					
Level	Conc.	Calc Conc.	CPS	Det.	RSD
1	0	0	1 P		86.6
2	10	9.89	13676 P		0.6
3	50	49.984	69112 P		0.9
4	100	100.019	138295 P		0.4

# NEW 7900 Ultra High Matrix Introduction (UHMI)


- Increased dilution range to **x100** – even higher matrix capability
- **Less matrix loading** to interface, so better long-term stability
- Maintain **high carrier gas flow** through spray chamber, so faster gas replacement and washout

	7700 HMI	New UHMI
HMI-4 (HMI-L)	0.6 L/min	0.8 L/min
HMI-8 (HMI-M)	0.35	0.68
HMI-25 (HMI-H)	<b>0.23</b>	<b>0.5</b>
<b>HMI-50</b>	N.A.	0.4
<b>HMI-100</b>	N.A.	0.33



# Matrix Tolerance Vastly Superior to Any Other ICP-MS

- Analysis of Brine (**25% NaCl**) with **UHMI (HMI-50)** & **new Humidifier**





- 9 Be [ no ]
- 27 Al [ no ]
- 51 V [ no ]
- 51 V [ He ]
- 52 Cr [ no ]
- 52 Cr [ He ]
- 55 Mn [ no ]
- 55 Mn [ He ]
- 56 Fe [ He ]
- 59 Co [ no ]
- 59 Co [ He ]
- 60 Ni [ no ]
- 60 Ni [ He ]
- 63 Cu [ He ]
- 66 Zn [ no ]
- 66 Zn [ He ]
- 68 Zn [ no ]
- 68 Zn [ He ]
- 75 As [ He ]
- 78 Se [ He ]
- 107 Ag [ no ]
- 107 Ag [ He ]
- 111 Cd [ no ]
- 111 Cd [ He ]
- 205 Tl [ no ]
- 205 Tl [ He ]
- 208 Pb [ no ]
- 208 Pb [ He ]
- 238 U [ no ]
- 238 U [ He ]

	OES DL (µg/L)	7900 DL (µg/L)
Ag	0.7	<b>0.07</b>
Al	0.5	0.6
Be	0.07	<b>0.03</b>
Cd	0.4	<b>0.16</b>
Co	1	<b>0.04</b>
Cr	1	<b>0.05</b>
Fe	1.7	<b>0.33</b>
Mn	0.2	<b>0.12</b>
Ni	2.5	<b>0.12</b>
Pb	6	<b>0.73</b>
Se	10	<b>1.2</b>
Tl	9	<b>0.22</b>
V	1	<b>0.09</b>
Zn	0.5	<b>0.39</b>

- Better DL than OES
- No maintenance required for cones and lenses after 8 hours continuous introduction.

# Agilent Ar gas humidifier

- Hollow fiber membrane tubing
- 2 channels in one body
  - For Carrier gas and Dilution gas
- No delay on gas flow change by much lower volume
  - Easier to tune (incl. Auto tune)



# NEW 7900 Orthogonal Detector System

Improved signal to noise and dynamic range

## Higher sensitivity

- High sensitivity EM (increase secondary electron generation by higher voltage at the 1st dynode)

## Lower background


- Off-axis from Q-pole to Detector

## Improved S/N


- Reduced noise on pulse signal
- New advanced discriminator system – separate noise by pulse width and pulse height

## Extended upper limit of detector range

- Analog range increased from 3 Gcps → 10 Gcps
- Up to 11 orders dynamic range - 0.1 cps to 10 Gcps




# Time Domain Discriminator



# Far Wider Measurement Range Than Any Other ICP-MS


11 orders - low and high level calibrations in a single run

**Cd (1ppt - 1ppb) and Na (100ppb - 10,000ppm (1%))** in the same run



Analyte Information (111 Cd [No Gas])							
	Current Sample	Calc Conc.	CPS	Ratio	Det.	Conc. RSD	
	1ppbBeAsCdHgPb	1000.006	50630.23		P	0.4	
Calibration							
Level	Rjct	Conc.	Calc Conc.	CPS	Ratio	Det.	RSD
1	<input type="checkbox"/>	0.000	0.000	3.00		P	47.1
2	<input type="checkbox"/>	1.000	0.960	51.60		P	14.9
3	<input type="checkbox"/>	10.000	10.078	513.21		P	3.8
4	<input type="checkbox"/>	50.000	49.851	2526.83		P	2.9
5	<input type="checkbox"/>	100.000	100.008	5066.11		P	2.2
6	<input type="checkbox"/>	1000.00	1000.006	50630.23		P	0.4
7	<input type="checkbox"/>						
8	<input type="checkbox"/>						
9	<input type="checkbox"/>						
10	<input type="checkbox"/>						
11	<input type="checkbox"/>						
12	<input type="checkbox"/>						
13	<input type="checkbox"/>						
14	<input type="checkbox"/>						

Both calibrations are linear.  
 Total concentration range covered from Cd blank (BEC of <0.1ppt) to Na top standard (1%) is 11 orders




Analyte Information (23 Na [HMI He] ISTD: 115 In [HMI He])							
	Current Sample	Calc Conc.	CPS	Ratio	Det.	Conc. RSD	
	rins	55.072	3281161.01	52.9048	A	31.7	
Calibration							
Level	Rjct	Conc.	Calc Conc.	CPS	Ratio	Det.	RSD
1	<input type="checkbox"/>	0.000	0.000	1944.19	0.0477	P	5.8
2	<input type="checkbox"/>	0.100	0.109	6393.55	0.1528	P	2.9
3	<input type="checkbox"/>	1.000	1.049	43220.26	1.0541	P	1.3
4	<input type="checkbox"/>	10.000	10.192	421823.28	9.8296	P	3.2
5	<input type="checkbox"/>	100.000	99.124	4004216.88	95.1848	A	1.7
6	<input type="checkbox"/>	1000.00	1089.335	40455458.40	1,045.570	A	1.5
7	<input type="checkbox"/>	10000.0	9991.075	328778051.1	9,589.288	A	3.7
8	<input type="checkbox"/>						
9	<input type="checkbox"/>						
10	<input type="checkbox"/>						
11	<input type="checkbox"/>						
12	<input type="checkbox"/>						
13	<input type="checkbox"/>						
14	<input type="checkbox"/>						

Concentration range (11 orders) and upper measurement limit (>1%) are at least 10x better than any other ICP-MS

# Faster TRA mode for Transient Signal Analysis Required for Single Nanoparticle analysis

Minimum dwell time is shortened from **3 msec.** on the 7700 Series to **0.1 msec** to allow faster sampling of transient signals.



**30 nm Au NP**

Acq Method | Data Analysis Method | Sample List

Acq Parameters | PeriPump/ISIS | Tune

**Acq Mode**

TRA

**Acq Option**


Tune Mode		#1: No Gas
Mass	Element Name	Monitor
197	Au	<input checked="" type="checkbox"/>

IntegTime /Mass [sec]
0.0001



# Faster TRA mode for Transient Signal Analysis


Minimum dwell time is shortened from **3 msec.** on the 7700 Series to **0.1 msec** to allow faster sampling of transient signals.




Time Resolved Analysis signal for single 30nm Au nanoparticle (dwell time 0.1ms)

# ORS<sup>4</sup> - Rapid ORS (Fast cell gas switching)

**7700**



**7900**




**7700**

Tune mode	No gas	He	HEHe	H2
Stabilization time	0 sec	10 sec	10 sec	15sec


**7900**

Tune mode	No gas	H2	He	HEHe
Stabilization time	0 sec	5 sec	5 sec	5sec

# Fast Cell Gas Switching – ORS<sup>4</sup>



# New Integrated Sample Introduction System (ISIS 3)





## New features in ISIS 3

- **Close-coupled valve** – very short tube length so minimal stabilization/rinse delay
- Easier tubing setup by **color coding**
- **Piston pump** for faster sample uptake
- **3-way valve** to switch between on-line ISTD or tune solution
  - ISIS is now compatible with Startup auto-optimization functions and full autotune

# Tune/ISTD valve for ISIS 3

- Support Start up, Autotune in batch
- Switch Tuning solution and ISTD solution automatically



# Improved Maintenance and Usability

Designed to make maintenance simpler and faster



U-series connection for MicroMist

Tighter connection S/C and Endcap

Peristaltic Pump redesigned:

- Easier to connect tubing, no drip onto Peri-pump roller

Torch box cover now uses clip (not thumb screw)




# Analysis of Major and Trace elements in Brine

## Sample

- Saturated salt water solution (~23% NaCl)
  - Typically analyzed by other than ICP-MS such as ICP-OES

## Experiment

- Agilent 7900 ICP-MS with UHMI
- Ar gas humidifier
- HMI-25



# Instrument Parameters

Mode	No gas	H <sub>2</sub>	He	High Energy He
Stabilization, sec	0	5	5	5
RF Power, W	<b>Preset Plasma</b> <b><u>HMI-25</u></b>			1600
Sampling Depth, mm				10
Carrier gas, L/min				0.45
Dilution gas, L/min				0.36
Ext 2, V	-235	-225	-220	-195
Omega Bias, V	-110	-120	-120	-100
Omega Lens, V	10.6	9.8	9.6	8.6
Deflect, V	10.8	0.6	1.0	-70
ORS gas, mL/min	-	6	5	10


**Quick stabilization**

**Ion lens**  
**Auto tuned**



# Recovery test

Very good recoveries are obtained over a wide concentration range from 20ppb to 250ppm.




Spiked conc.:

K	250 ppm	Sr	0.4 ppm
S	62.5 ppm	Al, Ti	0.2 ppm
Br	50 ppm	Others	20 ppb
Ca	40 ppm		

# Long term stability (5 hours)


ISTD correction works very well, and shows good signal stability over 5 hours even such a high TDS sample.



No need to clean lenses after running brine

# Analysis of various type of table salts

Covers wide measurement range from <math><0.01\text{ ppb}</math> to >1,000 ppm in Brine.



Salt samples #1-15:

Sea or Rock salt from Japan, Mexico, Germany, S.Africa, USA, Pakistan, Mongol

# Determination of Trace Elements in Steel

## Sample

- 17 Steel CRM from Japan Iron and Steel Federation (JISF)

## Sample prep.

- Dissolving CRM with HCl/HNO<sub>3</sub>
- Final concentration of Fe in solution: **0.5%** (JISF method)

## Instrument

- Agilent 7900 ICP-MS with UHMI
- HMI-25

# Analytical Conditions


## Preset Plasma HMI-25

	No gas	H <sub>2</sub>	He	HEHe
RF Power	1600 W			
Sampling Depth	10 mm			
Carrier Gas	0.45 L/min			
Dilution Gas	0.54 L/min			
Ext 1	-5	-5.3	-2.5	-7
Ext 2	-180	-195	-200	-200
Omega B	-105	-120	-115	-115
Omega L	9.6	10.9	10.7	7.2
Cell Ent	-36	-34	-26	-150
Cell Exit	-60	-70	-70	-150
Deflect	11.8	-1.4	-1.6	-75
Plate B	-45	-60	-60	-150
ORS Gas		5 mL/min	4 mL/min	10 mL/min
Oct P Bias	-8	-20	-20	-100

## Ion lens Auto tuned



## Cell parameters Preset

# ISTD Signal Stability (0.5% Steel, HMI-25)





Signal drift is **less than 20%** for most of ISTD elements even after **6 hours** even running **0.5% dissolved solution!**

# Relationship between Counts (corrected by ISTD) and Certified Value



# Relationship between Counts (corrected by ISTD) and Certified Value





# High matrix samples (Soil) analysis with ISIS 3

- Maximizing productivity



# ISIS 3 Program on MassHunter 4.1

Batch - 20131129\_ISIS(Soil)Run10.b

Save Batch Add to Queue Validate Method

Acq Method Data Analysis Method Sample List

Acq Parameters PeriPump/ISIS Tune

Sample Introduction: ISIS Discrete Sampling Tune Vial:

Sample Acquisition						
	Time [sec]	Speed [rps] Nebulizer Pump	Speed [%] ISIS Uptake Pump	Vial#	Valve	
<b>Common</b>						
Carrier Speed						
<b>Pre Run</b>						
Sample Load	6	Tune Parameter	50	Sample	Load	
Stabilize	7	Tune Parameter	5	Rinse Port	Inject	
<b>Acquisition</b>						
Probe Rinse (Sample)	23	Tune Parameter	5	Rinse Port	Inject	
Probe Rinse (Std)	23	Tune Parameter		Rinse Port	Inject	
Rinse 1		Tune Parameter	80		Inject	
Probe Rinse 1	5	Tune Parameter		Rinse Port	Inject	
Rinse 2		Tune Parameter	5		Inject	
Probe Rinse 2	10	Tune Parameter		Rinse Port	Inject	
Rinse 3		Tune Parameter			Inject	
Probe Rinse 3		Tune Parameter		Rinse Port	Inject	
<b>PostRun</b>						
Opt Loop Probe Wash	7	Tune Parameter	50	Rinse Port	Load	

Intelligent Rinse

**ISIS 3 pump program** enables to set optimized ISIS pump program to reduce carry over **during rinse steps.**



- Slow down the pump speed to fill the rinse solution in the rinse port.
- Wash the outside and inside of autosampler probe.

# Cycle time

## 83 sec. / sample

- Water, Soil, Sediment CRM
- 26 elements (ppb - ppm)
- **>20% faster** than 7700 ISIS (App. Note 5990-5437EN)

078SMPL.d	2013-11-11 00:23:51	Sample	50x Soil B
079SMPL.d	2013-11-11 00:25:12	Sample	1640a
080SMPL.d	2013-11-11 00:26:35	Sample	10x 1640a
081SMPL.d	2013-11-11 00:27:57	Sample	10x SW
082SMPL.d	2013-11-11 00:29:21	Sample	10x SW MS
083SMPL.d	2013-11-11 00:30:46	Sample	10x SW MSD
084_QC1.d	2013-11-11 00:32:10	QC1	LLCCV
085_QC2.d	2013-11-11 00:33:32	QC2	CCV
086_QC3.d	2013-11-11 00:34:55	QC3	CCB
087SMPL.d	2013-11-11 00:36:18	Sample	10x RS-A
088SMPL.d	2013-11-11 00:37:38	Sample	50x RS-A
089SMPL.d	2013-11-11 00:38:59	Sample	10x RS-B
090SMPL.d	2013-11-11 00:40:20	Sample	50x RS-B
091SMPL.d	2013-11-11 00:41:43	Sample	10x ES
092SMPL.d	2013-11-11 00:43:04	Sample	50x ES
Sample			10x Soil-A
Sample			50x Soil A
Sample			10x Soil B
Sample			50x Soil B



# Method Detection Limit

	MDL (ppb)
9B	0.07
23Na	6.1
24Mg	1.3
27Al	1.3
39K	7.8
44Ca	12
51V	0.03
52Cr	0.04
55Mn	0.04
56Fe	0.74
59Co	0.02
60Ni	0.05
63Cu	0.02

	MDL (ppb)
66Zn	0.06
75As	0.02
78Se	0.25
98Mo	0.02
107Ag	0.01
111Cd	0.04
121Sb	0.02
137Ba	0.05
201Hg	0.01
205Tl	0.01
208Pb	0.03
232Th	0.02
238U	0.01


# Quant Results



## 10x Diluted Sediment and Soil CRMs (ng/g)

	RS-A		RS-B		ES		Soil-A		Soil-B	
	measured	certified	measured	certified	measured	certified	measured	certified	measured	certified
Be	0.006		0.013		1.94	2.00	0.011		0.010	
Na	5297	5000	5012	5000	19770	20000	7246	7000	10090	10000
Mg	6897	7000	11960	12000	10030	10000	6960	7000	8189	8000
Al	25940	25000	59500	60000	68820	70000	50480	50000	68010	70000
K	15001	15000	20030	20000	14750	15000	19510	20000	20230	21000
Ca	29950	30000	29340	30000	7927	8000	34320	35000	11960	12500
V	27.5	25.0	105.4	100	103.6	100	10.0	10.0	80.2	80.0
Cr	30810	30000	1521	1500	81.4	80.0	0.192		39.9	40.0
Mn	876	800	612	600	407	400	10.4	10.0	9665	10000
Fe	127800	120000	40280	40000	35730	35000	19550	20000	34610	35000
Co	12.1	10.0	15.7	15.0	10.7	10.0	0.317		10.5	10.0
Ni	52.8	50	51.1	50	30.4	30	29.4	30	20.7	20
Cu	102.6	100	100.0	100	20.1	20.0	28.7	30.0	300	300
Zn	1531	1500	497	500	151	150	96.8	100	6798	7000
As	59.4	60	19.5	20	9.90	10.0	19.8	20	597	600
Se	2.41	2.0	1.13	1.0	4.89	5.0	0.925	1.0	0.089	
Cd	11.1	10	3.37	3.0	0.099		0.347		19.4	20.0
Sb	71.6	50	4.38	4.0	0.548		3.12	3.0	39.4	40.0
Ba	54.9	50	421	400	1.30		493	500	693	700
Tl	1.01	1.0	1.18	1.0	0.014		0.022		0.179	
Pb	714	700	202	200	29.8	30	39.8	40	5793	6000
Th	1.94	2.0	10.2	10.0	10.2	10.0	10.1	10.0	10.2	10.0
U	0.94	1.0	3.01	3.0	0.001		0.979	1.0	25.4	25

# CCV Stability



Na, Mg, K, Ca, Fe: 5 µg/g  
 Other elements: 50 ng/g

# Interface

### Sampler





### Skimmer



### Tip of Skimmer



**Before**  
conditioning  
and  
sequence



**After**  
conditioning  
and  
sequence

# Improving Installation & Familiarization

## Installation


- Intuitive Installation Checkout
- Site Prep tool for Software Installation

## Familiarization

- Familiarization Tutorials
  - Familiarization Guide
  - **Familiarization Video (over 20 video clips for operation)**
  - Familiarization Slide Set

Plus:

**Remote Advisor** Now available for ICP-MS!





# Summary – New Agilent 7900 ICP-MS


## Better customer experience

### **Better Analytical performance experience**


- Ultra high matrix tolerance (UHMI)
- Superior sensitivity and lower background noise (ODS)
- Wider dynamic range (ODS)
- New Productivity Option (ISIS 3)
- Ultra fast scan speed for Single Nanoparticle analysis

### **Better Support experience**

- Familiarization Tutorials/Videos
- Remote Advisor support



# New Agilent 7900 ICP-MS



**New Plasma interface & Vacuum design**  
Enables higher **Sensitivity**



**New Orthogonal Detector system**

- Higher **Sensitivity**
- Lower **Back ground**
- Extended Upper limit



**New UHMI**  
Enables analysis of 25% brine

**ORS<sup>4</sup>**

3 sec. cell gas switching

**New MassHunter 4.1**

- Improved User Experience
- Method automation
- SPSP common reporting engine

**New Processor card**  
**New Firmware**

0.1msec dwell time for fast transient and single nanoparticle analysis

**PLUS**

- ✓ New Humidifier
- ✓ Improved Usability
- ✓ 2 Days Installation
- ✓ Operating Video tool
- ✓ Remote Advisor
- ✓ Less Maintenance

**New Styling**

Consistent across SPSP products

**New Productivity option (ISIS 3)**

Faster and Fully automated from Start-up

